

PIXEL TECHNOLOGY

The Ultimate Hi-Tech Eyes
and Lips Corrector



GP
General Project
ADVANCED BIOMEDICAL TECHNOLOGIES

CRYO
FACE RF

www.generalproject.com

PIXEL

TECHNOLOGY

General Project is proud to present its latest endeavour in the field of Research and Innovation, the newest Facial Rejuvenation device, **CryoFace RF**, exploiting the ultimate Pixel Technology energy emission.

CryoFace RF takes advantage of the strength and efficiency of two combined technologies - Hyper Pulsed RF and Cooling - to fight the most annoying signs of facial ageing, as crow's feet and lips barcode, with a treatment that effectively smoothes wrinkles and reduces skin laxity.

Energy Booster handpiece



The innovative **Pixel Technology** is driven by the special **Energy Booster handpiece** that, thanks to its specific size and shape, focuses precisely on the periocular and perilabial areas, delivering in depth high quantity of energy duly counterbalanced by cooling. This special handpiece can also treat directly the whole lips surface, thus acting also as a volume enhancer.

One further RF standard handpiece is available to perform complete treatments on large areas, as face, neck and décolletage, to immediately improve skin texture, tone and clarity.

TECHNICAL FEATURES

Power Supply:	220 ÷ 240 Volt AC (100 ÷ 130V) 50/60 Hz single phase
RF Frequency:	1 MHz (±20 %)
Cooling:	down to 2°C

GP
General Project
ADVANCED BIOMEDICAL TECHNOLOGIES

General Project s.r.l.

Via della Gora 15/19 - 50025 Montespertoli (Firenze) Italy - Tel. +39 0571 675076 - Fax +39 0571 675077
www.generalproject.com - info@generalproject.com